Temperierbad Additiv

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REF 200 001: 100 ml

Temperierbad Additiv

Concentrated microbiocidal additive for water soluble heating bath

Application and compatibility

Concentrated microbiocidal additive for clean keeping of water soluble heating bath.

This Additive inhibits the growth of bacteria, algae, fungi and inhibits the precipitation of inorganic salts.

This product is validated for the use with water soluble heating bathes.

This Additive is **not compatible** with **aluminium & acryl-glas** materials.

The application has to be performed in reference to this package insert released by **invicon**, otherwise no guarantee is valid.

This product is intended to be used only by laboratory personnel according to the guidelines for laboratories (GLP).

This package insert contains important information. Please read it carefully and keep it.

Handling

Add 5 ml of the additive to 1 l water bath solution and mix. Renew the water bath solution and the additive every 4 weeks. In the case of blue color of solution faded, please renew earlier.

Contamination-Indicator

Bacterial contamination is indicated by decoloration of the solution.

Storage and stability

Reagents are stable up to the declared expiration date when stored at 18...25 °C.

Reagent components, active components

Quarternary ammonium compounds, complexing agents, indicator, surfactants
Free of aldehydes and azides

Precautions, warnings

H302 Harmful if swallowed.

H314 Causes severe skin burns and eye damage.

H400 Very toxic to aquatic life.

- Store the reagents protected from children.
- This product is to be used only by special trained laboratory personnel. Exercise the usual precautions required for handling all laboratory reagents.
- Do not swallow, do not pipette with the mouth, avoid contamination and contact with the skin. Do not inhale.

Disposal instructions

Additive is biological degradable. Packagings and completely empty and rinsed bottles or vials are to be disposed according to the national guidelines.