

# Seronorm™ CRP Liquid

Straight from the fridge into the instrument

≥ 95%  
human  
serum

## Analytes:

CRP

Refer to the lot specific product documentation of currently available lots for specific analyte and stability claims. For details concerning intended use, please refer to the Instructions for Use (IFU). The documents are available on our website [www.sero.no](http://www.sero.no) or upon request.



## Main features

Contains more than 95% human serum

Liquid material stable at 2–8°C

Three clinically relevant levels

1-month stability in opened vials (2–8°C)

Shelf life of 3 years

Independent QC material

## Benefits

*Behaves similar to a patient sample (high commutability), allowing highly trustworthy results, especially during reagent lot changes.*

*Straight from the fridge to the instrument. Ready to use; no thawing or reconstitution needed. Eliminates potential handling errors.*

*Ensures patient safety with a broader quality control range.*

*Reduces unnecessary waste for small-volume users.*

*Longer shelf life that makes it possible to reduce the frequency of lot change: resource and time-related consequences.*

*Provides an objective assessment of the entire analytical system.*

## Stability

Shelf life	Open vial stability	
	2–8 °C	15–25 °C
2–8 °C	2–8 °C	15–25 °C
3 years	31 days	7 days

## Order Information:

Art. no	Product	Size
213005	Seronorm™ CRP Liquid L-1	12 x 1 mL
213105	Seronorm™ CRP Liquid L-2	12 x 1 mL
213205	Seronorm™ CRP Liquid L-3	12 x 1 mL