



## DHEA Plus

UTAK recognizes the need for laboratories to have accurate and reliable test results in the steroid hormone field so we have developed the highest quality Human Serum Quality Controls for steroid hormones.

The steroid or sex hormones may be the first that come to mind for most people. However, just as significant are the steroid hormones produced by the adrenal glands. This group of hormones supports a broad range of essential physiological functions.

According to Women's Health Connection, "Optimal hormone balance is critical to our physical and emotional well-being, as well as long-term protection against age-related concerns."

UTAK, beyond measure... one test result at a time

### 100% REAL Serum

- Reduce Matrix Effects
- Increase Laboratory Efficiency

### Open Architecture 3rd Party Control

- Increase Laboratory Flexibility
- Use on Any Testing Platform



# DHEA Plus Serum Control



## Specifications:

Analyte	Low	High
	REF 51410	REF 51411
	Reference Value	Reference Value
Corticosterone	3 ng/mL	8 ng/mL
Cortisol	60 ng/mL	200 ng/mL
Cortisone	15 ng/mL	40 ng/mL
Dehydroepiandrosterone (DHEA)	3 ng/mL	8 ng/mL
Progesterone	7.5 ng/mL	20 ng/mL
17 $\alpha$ -Hydroxyprogesterone	3 ng/mL	8 ng/mL
Testosterone	3 ng/mL	8 ng/mL
Androstenedione	1.5 ng/mL	4 ng/mL
11-Deoxycortisol	100 ng/dL	300 ng/dL

Form: Frozen

Storage: -10° C (14° F)

Configuration: 51410 Low 5x3 mL  
51411 High 5x3 mL

Shelf Life: 12 Months from Manufacture Date

Reconstituted  
Stability: 25 Days

For more information, please contact your distributor.

Service und Vertrieb Deutschland:



25020 Avenue Tibbitts  
Valencia, CA 91355



Scan Here to  
Visit [UTAK.com](http://UTAK.com)