



## **Steroids Panel**

UTAK recognizes the need for laboratories to have accurate and reliable test results in steroid hormone field, so we have added to our portfolio, a new Quality Control for steroid hormones. These expanded level controls are manufactured in our new proprietary Specialty Matrix SMx<sup>™</sup> Serum.

The steroid or sex hormones may be the first that come to mind for most people. However, just as significant are the steroid hormones produced by the adrenal glands. This group of hormones supports a broad range of essential physiological functions.

According to Women's Health Connection, "Hormone testing is a valuable tool for obtaining optimal hormone balance, as well as for maintaining that balance once it is achieved."


UTAK, beyond measure...one test result at a time.

## SMx™ Serum

- Reduce Matrix Effects
- Increase Laboratory Efficiency

## **Open Architecture 3rd Party Control**

- Increase Laboratory Flexibility
- Use on Any Testing Platform





Specifications:		Negative	Level 1		Level 2	
		<b>REF</b> 51400	REE 51401	[	REE 51402	
Analyte		Reference Value	Reference Val	ue Re	ference Valu	Ie
Cortisol		0 ng/mL	25 ng/mL		50 ng/mL	
Cortisone		0 ng/mL	6 ng/mL		12 ng/mL	
Corticosterone		0 ng/mL	1 ng/mL		2 ng/mL	
DHEA		0 ng/mL	1 ng/mL		2 ng/mL	
Progesterone		0 ng/mL	2.5 ng/mL		5 ng/mL	
17α-Hydroxyprogesterone		0 ng/mL	1 ng/mL		2 ng/mL	
Androstenedione		0 ng/mL	.5 ng/mL		1 ng/mL	
11-Deoxycortisol		0 ng/mL	40 ng/dL		80 ng/dL	
Testosterone		0 ng/mL	1 ng/mL		2 ng/mL	
		Level 3	Level 4			
		<b>REF</b> 51403	<b>REF</b> 51404			
Analyte		Reference Value	Reference Val	ue		
Cortisol		125 ng/mL	250 ng/mL			
Cortisone		30 ng/mL	60 ng/mL			
Corticosterone		5 ng/mL	10 ng/mL			
DHEA		5 ng/mL	10 ng/mL			
Progesterone		12.5 ng/mL	25 ng/mL			
17α-Hydroxyprogesterone		5 ng/mL	10 ng/mL			
Androstenedione		2.5 ng/mL	5 ng/mL			
11-Deoxycortisol		200 ng/dL	400 ng/dL			
Testosterone		5 ng/mL	10 ng/mL			
Form:	Frozen			Storag	e: -10	°C (14°F)
Configuration:		egative 5x3 m		Shelf L	_ife: 12 I	Months from Manufactu
		evel 1 5x3 m				
		evel 2 5x3 m		Thawe		<b>D</b>
		evel 3 5x3 m		Stabili	ty: 25 l	Days
	51404 Le	evel 4 5x3 m	۱L			

For more information, please contact your distributor.







Scan Here to Visit UTAK.com