



## Amphetamines Panel

Amphetamine abuse has skyrocketed across the globe with the popularity of MDMA or “Ecstasy” and the highly addictive and toxic “Crystal Meth” or methamphetamine. Abuse of Amphetamines can cause addiction, neurotoxicity, depression, and death.

We make our Amphetamine Quality Controls with 100% Real Human Urine to give laboratorians the assurance of the most accurate test results.

UTAK, beyond measure...one test result at a time.

### **100% Real Human Urine**

- Reduce Matrix Effects
- Increase Laboratory Efficiency

### **Open Architecture 3rd Party Control**

- Increase Laboratory Flexibility
- Use on Any Testing Platform



# Amphetamines Panel Urine Control



## Specifications:

Analyte	Level 1	Level 2
	REF 41400	REF 41401
	Reference Value	Reference Value
Amphetamine	350 ng/mL	700 ng/mL
MDA	350 ng/mL	700 ng/mL
Methamphetamine	350 ng/mL	700 ng/mL
MDMA	350 ng/mL	700 ng/mL
MDEA	350 ng/mL	700 ng/mL
Ephedrine	350 ng/mL	700 ng/mL
Methylephedrine	350 ng/mL	700 ng/mL
Phendimetrazine	350 ng/mL	700 ng/mL
Phenmetrazine	350 ng/mL	700 ng/mL
Phentermine	350 ng/mL	700 ng/mL
Phenylpropanolamine	350 ng/mL	700 ng/mL
Pseudoephedrine	350 ng/mL	700 ng/mL
Selegiline	350 ng/mL	700 ng/mL

**Form:** Lyophilized

**Storage:** 2-8°C (35-46°F)

**Configuration:** 41400 Level 1 5x5 mL  
41401 Level 2 5x5 mL

**Shelf Life:** 30 Months from Manufacture Date

**Reconstituted  
Stability:** 25 Days

For more information, please contact your distributor.

Service und Vertrieb Deutschland:

**invicon**  
diagnostic concepts GmbH  
Agnes-Pockels-Bogen 1 • 80992 München • Germany  
+49 89 319047-0 • info@invicon.de • www.invicon.de

**UTAK**  
createCONTROL

25020 Avenue Tibbitts  
Valencia, CA 91355



Scan Here to  
Visit [UTAK.com](http://UTAK.com)