



## Miscellaneous Panel

Behind each test result is a real life human being with a story to tell. Drug test results change lives, which is why UTAK Miscellaneous Panel Quality Controls are made with only 100% Real Human Urine to give laboratorians the most accurate and reliable results.

UTAK's ongoing commitment to the manufacture of 100% Real Human Matrix Controls is a commitment to the individual lives of those relying on accurate and reliable test results.

UTAK, beyond measure... one test result at a time

### **100% Real Human Urine**

- Reduce Matrix Effects
- Increase Laboratory Efficiency

### **Open Architecture 3rd Party**

- Increase Laboratory Flexibility
- Use on Any Testing Platform



# Miscellaneous Panel Urine Control



## Specifications:

Analyte	Level 1	Level 2
	REF 41430	REF 41431
	Reference Value	Reference Value
Carisoprodol	500 ng/mL	1,500 ng/mL
Cocaine	100 ng/mL	300 ng/mL
Benzoyllecognine	100 ng/mL	300 ng/mL
Cocaethylene	100 ng/mL	300 ng/mL
Cyclobenzaprine	500 ng/mL	1,500 ng/mL
Meprobamate	500 ng/mL	1,500 ng/mL
Methaqualone	100 ng/mL	300 ng/mL
Methylphenidate	500 ng/mL	1,500 ng/mL
Phencyclidine (PCP)	100 ng/mL	300 ng/mL
Ritalinic Acid	500 ng/mL	1,500 ng/mL
Zolpidem	200 ng/mL	600 ng/mL
(l)-9-Carboxy-11-nor- $\Delta^9$ -THC	100 ng/mL	300 ng/mL

**Form:** Lyophilized

**Storage:** 2-8°C (35-46°F)

**Configuration:** 41430 Level 1 5x5 mL  
41431 Level 2 5x5 mL

**Shelf Life:** 30 Months from Manufacture Date

**Reconstituted  
Stability:** 25 Days

For more information, please contact your distributor.

Service und Vertrieb Deutschland:



25020 Avenue Tibbitts  
Valencia, CA 91355



Scan Here to  
Visit [UTAK.com](http://UTAK.com)