

# ETHANOL



It's this straightforward, UTAK 100% human matrix controls are the standard against which all others are judged.

## **UTAK** ETHANOL CONTROLS:

- ✓ OFFER THE MOST RELIABLE, CONSISTENT RESULTS
- ✓ DELIVER THE BEST VALUE IN TIME AND MONEY SAVED
- ✓ ARE BACKED BY UNMATCHED CUSTOMER SERVICE



**UTAK**<sup>®</sup>  
LABORATORIES, INC.  
*the leader in 100% human matrix controls*

## BI-LEVEL ETHANOL Whole Blood Toxicology Controls

### ETHANOL 50 mg/dL

Analyte	Target Value
Ethanol	<b>50 mg/dL</b> (10.8 mmol/L)

Form: Frozen

Shelf Life: 12 Months

Volume: 5x3 mL

### ETHANOL 80 mg/dL

Analyte	Target Value
Ethanol	<b>80 mg/dL</b> (17.4 mmol/L)

Form: Frozen

Shelf Life: 12 Months

Volume: 5x3 mL



**UTAK**<sup>®</sup>  
LABORATORIES, INC.  
*the leader in 100% human matrix controls*