



## Benzodiazepines Plus

Benzodiazepines are often prescribed as depressant medications. Many different types of benzodiazepine medications exist to treat both psychological and physical conditions.

Consistent and accurate monitoring for compliance and to identify misuse and abuse of these medications is essential toward making appropriate treatment decisions.

UTAK's Handcrafted Benzodiazepines Plus Quality Control affords accuracy to the discipline of analytical measurement.

UTAK, beyond measure... one test result at a time

### **Open Architecture 3rd Party Control**

- Increase Laboratory Flexibility
- Use on Any Testing Platform

### **100% Real Human Urine and Whole Blood**

- Reduce Matrix Effects
- Increase Laboratory Efficiency



# BENZODIAZEPINES PLUS

## Urine Controls & Whole Blood Controls



### Specifications:

Analyte	URINE 100	URINE 400	WB 100
	<small>[REF]</small> 12090 Reference Value	<small>[REF]</small> 12091 Reference Value	<small>[REF]</small> 12100 Reference Value
Alprazolam	100 ng/mL	400 ng/mL	100 ng/mL
α-Hydroxyalprazolam	100 ng/mL	400 ng/mL	100 ng/mL
Chlordiazepoxide	100 ng/mL	400 ng/mL	100 ng/mL
Norchlordiazepoxide	100 ng/mL	400 ng/mL	n/a
Clobazam	100 ng/mL	400 ng/mL	100 ng/mL
Clonazepam	100 ng/mL	400 ng/mL	100 ng/mL
7-Aminoclonazepam	100 ng/mL	400 ng/mL	100 ng/mL
Diazepam	100 ng/mL	400 ng/mL	100 ng/mL
Nordiazepam	100 ng/mL	400 ng/mL	100 ng/mL
Estazolam	100 ng/mL	400 ng/mL	100 ng/mL
Flunitrazepam	100 ng/mL	400 ng/mL	100 ng/mL
7-Aminoflunitrazepam	100 ng/mL	400 ng/mL	100 ng/mL
Flurazepam	100 ng/mL	400 ng/mL	100 ng/mL
Desalkylflurazepam	100 ng/mL	400 ng/mL	100 ng/mL
Hydroxyethylflurazepam	100 ng/mL	400 ng/mL	100 ng/mL
Lorazepam	100 ng/mL	400 ng/mL	100 ng/mL
Midazolam	100 ng/mL	400 ng/mL	100 ng/mL
α-Hydroxymidazolam	100 ng/mL	400 ng/mL	n/a
Nitrazepam	100 ng/mL	400 ng/mL	100 ng/mL
7-Aminonitrazepam	100 ng/mL	400 ng/mL	100 ng/mL
Oxazepam	100 ng/mL	400 ng/mL	100 ng/mL
Temazepam	100 ng/mL	400 ng/mL	100 ng/mL
Triazolam	100 ng/mL	400 ng/mL	100 ng/mL
α-Hydroxytriazolam	100 ng/mL	400 ng/mL	100 ng/mL
Zolpidem	100 ng/mL	400 ng/mL	100 ng/mL

**Form:** Lyophilized

**Storage:** 2-8°C (35-46°F)

**Configuration:** 12090 Urine 100 5x5 mL  
12091 Urine 400 5x5 mL  
12100 WB 100 5x5 mL

**Shelf Life:** 30 Months from Manufacture Date

**Reconstituted Stability:** 25 Days

For more information, please contact your distributor.

Service und Vertrieb Deutschland:



25020 Avenue Tibbitts  
Valencia, CA 91355



Scan Here to  
Visit [UTAK.com](http://UTAK.com)